KP Digital Software Manual

1. Summary

KP Digital software is a software component of the K-Scan oral digital impression instrument and Smile Scan oral digital impression instrument produced by Guilin Kaiwen Peter Technology Co., Ltd. It runs independently in the Windows system to collect signals from the receiving host and input them into the software for model reconstruction and corresponding order management.

System Function Overview

KP Digital software consists of several modules, each of which performs specific functions to implement the case checking process and provide effective 3D data.

Functional Module	Sub modules/function points	Function point description
		users to register account information in the system and obtain account permissions
	log on	User enters account and password to log in
Registration and login	Remember password	You used to log in, so you don't need to enter it again the next time you log in
	Forgot password	Users reset their password in the system through the email bound during registration
	Import offline users	Users download the license file and import it into the system to import offline users
	Search Cases	Search based on patient name and case number and display search results
	Modify Case	Modifying Case Details
	New Case	Create a new case
Case	Delete Case	Delete selected case information
management	Upload case	Upload the selected case information and model files to the cloud
	Download Cases	Download model files from the cloud to the local location
	Copy Case	Copy the selected case
	Filter Cases	Filter cases based on their status and source
	Search patients	Search based on patient name and display search results
Patient Management	Modify Patient	Modifying Patient Details
	Create Patient	Create a new patient

The function introduction of each module of the software component is as follows:

1					
	Delete Patient	Delete selected patient information			
	Scan to obtain model data	Obtain current dental model data through scanning			
	Model optimization	Optimize the model data obtained through scanning to improve accuracy			
Model Scan	Complete Scan	Exit and save model data			
	Model processing	Can perform operations such as cutting, chamfering inspection, and drawing neck lines on the model			
	Model display processing	Can perform realistic display of models, prompt scanni quality, adjust model display size, take snapshots of mode and other operations			
	Model List	Display all models under the current case			
Model Preview	Model Display	Display the current 3D model screen			
	Angle adjustment	gle adjustment Switch the perspective of the current observation model			
	Search plugin	Search for plugins based on their names			
Plug-in	View plugin information	Click on the plugin to view the plugin details			
Manager	Installing plugins	Install the selected plugin to the local system			
	Uninstalling plugins	Uninstall the selected plugin			
	Commonly used	Software global language, folder name format for exporting cases, host can set storage path, backup path, and export path			
System settings	Help	Display manufacturer information, company website, and user manual			
	Software updates	Display the current software version and provide the function of manually checking the version, as well as the function of selecting offline upgrade packages for software updates			
Connection Assistant	Connection of mobile apps	Mobile APP scans QR code to connect to server			
User management	liser information	Display and edit the information of the logged in user, and can log out to the login interface			
Cloud management	Quickly jump to the cloud	Jump to the entrance of cloud web pages			

2 Operating environment

software environment

System software: Windows10 x64, Windows11 x64 **Network conditions** Network card: built-in network card

3 Software Function

3.1 Install

1

Follow the steps below to install the software.

Double click on the installation package to start the installation, as shown in the following figure

	Figure 1		
Select Setup Lang	guage	×	
Select the	language to use during the inst	allation.	
English		~	

Figure 2

Carefully read the license agreement, agree to it, click 'I agree to this agreement', click 'Next' to continue the installation, if you do not agree, click 'Cancel' to end the installation.

🛛 Setup - KP Digital version V1.0.0 —	×
License Agreement Please read the following important information before continuing.	ð
Please read the following License Agreement. You must accept the terms of this agreement before continuing with the installation.	
USE LICENCE FOR THE 《KP Digital》 SOFTWARE	1
THE PURCHASER, HEREINAFTER DESIGNATED AS THE « LICENSEE », EXPLICITLY ACCEPTS THE TERMS OF THE PRESENT LICENSE AND COMMITS HIMSELF THERETO IPSO FACTO BY THE MERE ACTION OF CLICKING ON THE «ACCEPTANCE» BUTTON, OR BY INSTALLING, COPYING OR USING IN ANY WAY WHATSOEVER THE "KP Digital" SOFTWARE, HEREINAFTER DESIGNATED AS THE «SOFTWARE».	
IF THE «LICENSEE» DOES NOT ACCEPT THE PRESENT AGREEMENT, HE MUST SELECT THE BUTTON INDICATING «NON-ACCEPTANCE» AND HE AGREES NOT TO INSTALL, COPY OR USE THE «SOFTWARE» IN ANY WAY WHATSOEVER.	
THE «LICENSEE» SHALL THEN HAVE TO CONTACT HIS RESELLER AS SOON AS POSSIBLE, WHICH WILL INFORM HIM ABOUT THE PROCEDURES FOR RETUNING THE «SOFTWARE».	
○ I <u>a</u> ccept the agreement	
I do not accept the agreement	
Next >	Cancel

Figure 3




Figure 5

Setup - KP Digital version V1.0.0 - X	
Select Additional Tasks Which additional tasks should be performed?	5
Select the additional tasks you would like Setup to perform while installing KP Digital, then click Next.	
Additional shortcuts:	
Create a desktop shortcut	
< <u>B</u> ack <u>N</u> ext > Cancel	
Figure 6	
Setup - KP Digital version V1.0.0 — 🗌 🔿	×
Ready to Install	
Setup is now ready to begin installing KP Digital on your computer.	À
Setup is now ready to begin installing KP Digital on your computer.	Z
	3
Click Install to continue with the installation, or click Back if you want to review or change any settings. Destination location: G: KP Digital Additional tasks: Additional shortcuts:	3
Click Install to continue with the installation, or click Back if you want to review or change any settings. Destination location: G: \KP Digital Additional tasks:	3
Click Install to continue with the installation, or click Back if you want to review or change any settings. Destination location: G: \KP Digital Additional tasks: Additional shortcuts:	Z
Click Install to continue with the installation, or click Back if you want to review or change any settings. Destination location: G: \KP Digital Additional tasks: Additional shortcuts:	
Click Install to continue with the installation, or click Back if you want to review or change any settings. Destination location: G: \KP Digital Additional tasks: Additional shortcuts:	A
Click Install to continue with the installation, or click Back if you want to review or change any settings. Destination location: G: \KP Digital Additional tasks: Additional shortcuts:	
Click Install to continue with the installation, or click Back if you want to review or change any settings. Destination location: G: \KP Digital Additional tasks: Additional shortcuts:	

Figure 7




Figure 8

Installation completed, click 'Finish'

Setup - KP Digital version V1.0.0	×
	Completing the KP Digital Setup Wizard Setup has finished installing KP Digital on your computer. The application may be launched by selecting the installed shortcuts. Click Finish to exit Setup.
	Einish

Figure 9

3.2 log on

After installation, there will be a software icon on the computer desktop, as shown in the figure. Double click the icon, and after the software starts, you will enter the software login interface, as shown in the figure.




Figure 11

3.2.1 Registration and login

If you have not registered, please click the registration button to go to the webpage to register your account, as shown in the figure.

lease select your account type	
Clinic	Laboratory
Clinic administrator With all software functions	Call Laboratory administrator With all software functions
Ø Order management	Order management
Organization member management	Organization member management
Partner management Register	Partner management Register
Clinic member	Laboratory member
Needs to be bound with the clinic, become a member	Needs to be bound with the technical institute,
of the clinic, and have the basic functions of the software	become a member of the technical institute, and have the basic functions of the software
Intraoral scanning	View order
View order Register	Ownload model file Register

Figure 12

Enter the username and password, click the login button to log in to the main interface, as shown in the figure. In the login interface, users can check "Remember Password" and do not need to enter the password again the next time they log in.

KP Digital V1.0.0	$\bigcirc English \lor \bigcirc \vspace{-1.5ex}$
	Welcome KP Digital
	Email 657587346@qq.com
KP P P	Password
	•••••
	Remember password Forgot password
	Sign in
	No account? Register

Figure 13

KP V1.0.0	Digital	R New Case	KP Digital Cloud Platform		G 🕀 🕯 🕬	O C
Case List	atient List Plugin List				Case No. , Patient Name	8
Case No.	Patient Name	Case Type	Case Status	Last Modified Time	Operation	
006	999	Repair	Waiting to scan	2023-08-19 16:08:45	🖸 🗊	
0041	bbbbbb酸噻噻嘎嘎噻唑剂则211	Repair	Waiting to deliver	2023-08-19 07:16:21	٩	
006	999	Repair	Waiting to scan	2023-08-18 14:17:58	🕒 🤠	
005	55	Repair	Waiting to scan	2023-08-18 11:41:10	🕒 🤠	
004	34	Repair	Waiting to scan	2023-08-18 11:16:48	()	
002	test	Repair	Accepted	2023-08-17 14:52:49	🖒 🖸 📋	
001	liangy	Repair	Accepted	2023-08-17 14:23:04		

7 in total Show per page 20 = (1/1)

Figure 14

If you forget your password, you can click "Forgot Password" on the login interface. The system will open the password retrieval interface on the webpage, as shown in the figure. After the user enters their email, the system will send an email. The user logs in to the email and opens the password reset link, where they can reset their password.

Reset password		
' Email		
Please input your er	nail	
[•] Captcha		
Please input image	captcha	ETW8C
	Send email	

Figure 15

If you have previously registered, you can download the license file in the personal center module of the KP Digital cloud platform, click on Import, select the license file to import, and then log in using the account, as shown in the figure




Figure 16

3.3 Case Management Module

3.3.1 Search Cases

Enter the case number or patient name in the input box, and then click the search icon. The system will display the corresponding search results, as shown in the figure

Case List Pati	ent List Plugin List				Case No. , Patient Name	8
Case No.	Patient Name	Case Type	Case Status	Last Modified Time	Operation	
006	999	Repair	Waiting to scan	2023-08-19 16:08:45	🕒 🔁	
0041	hhhhhh廣應環境環境開閉211	Repair	Waiting to deliver	2023-08-19 07:16:21	١	
006	999	Repair	Waiting to scan	2023-08-18 14:17:58	🖸 🔁	
005	55	Repair	Waiting to scan	2023-08-18 11:41:10	🖸 🗊	
004	34	Repair	Waiting to scan	2023-08-18 11:16:48	🔁 🔁	
002	test	Repair	Accepted	2023-08-17 14:52:49	🕹 🗈 📋	
001	liangy	Repair	Accepted	2023-08-17 14:23:04		



3.3.2 New Case

Click on the "New Case" icon on the homepage, as shown in the figure. Enter patient, dental position, and other information in the case interface, and then click Save, as shown in the figure \circ

1	KP Digital	New Case	KP Digital Cloud Platform		G 🔅 🛎 高仔	⊖¢
Case List	Patient List Plugin List				Case No. , Patient Name	8
Case No.	Patient Name	Case Type	Case Status	Last Modified Time	Operation	

0 in total Show per page 10 = (1/1)

Figure 18

KP Digital V1.0.0	•				5 6	💄 梁永	O C
< Case	Preview Send						Scan
Case Information		Tooth Selection					
ID	Create Time			Full Crow	n Missing Te	wth	Integ
009	2023-08-21 11:16:50		11 21	Venuor	Opposite 1	rech	
Case No.	Visit Type	12	22 23	Materials			
009	Initial diagnosis	14	24	Acrylic			
Patient Name*			Maxilla 25				
	L 🖸	16	26	Tooth Color			
Doctor*		17	27	A1			
梁永	~ 0	18	28	Preoperative Scar			
Case Type Repair O Plant	O Orthodontics O Others	(48)	38	No			
Remark	O Orthodontics O Others	47	37				
Special note: None.		H	H				
		46	Lower 36				
		45	35				
			33				
		×	2 41 31 32 33				
				Clear	Current	Clear	All
					Cancel	Browse	Save



3.3.3 Modify Case

Double-click the case to be modified in the case list to enter the case modification interface, as shown in the figure. Users can modify the corresponding information as needed. After modification, click the save button to save it $_{\circ}$

KP Digital V1.0.0 (D			FG 🕀 💄 🕯	^{₽\$} ⊖ C→
< Case	Preview Send				Scan
Case Information		Tooth Selection			
ID	Create Time		Full Cr	rown Missing Teeth	Intay
007	2023-08-19 16:08:45	(1) (21)	Vene	eer Opposite Teech	
Case No.	Visit Type	12 12 22			
006	Initial diagnosis		24 Materials		
Patient Name*		15 Maxilla	25 Acrylic		
999		16	26 Tooth Color		
Doctor*		17	27 A1		
梁永	- 0	18	28 Preoperative S	can	
Case Type			No		
Repair O Plant O	Orthodontics O Others	43	38		
Remark		47	37		
Special note: None.		40 43 44 41 41 41 41 41 41 41 41 41 41 41 41	35 35 34		
			Ch	lear Current	Clear All
				Cancel Brow	/se Save



3.3.4 Delete Case

Click the delete button in the case list, and the system will pop up a deletion prompt, as shown in the figure. After the user confirms the deletion, the selected case will be deleted and the case list will be refreshed.




Figure 21

3.3.5 Send Case

Click the send button in the case list to jump to the case sending interface, displaying the detailed information of the case, as shown in the figure. Select the technician to send, click the "send" button, and the data will be uploaded to the cloud.

KP Digital V1.0.0					F0 🕀 🔒	RA D C
Case Preview	Send					
Case Information	Tooth Selection	Teeth List				
D 002	11 21	Number	Type	Materials	Tooth Color	
Case No. 002	12 22	17	Full crown	Acrylic	A1	
Create Time 2023-08-17 14:46:15	14 24					
Doctor 梁永	15 Maxilla 25 16 26					
Patient test	17 27					
Case Type Repair	18 28					
Remark Special note: None.		Scanned Model			Send Settings	
	48 38				Target	
	47 37				Liangy_Lab	
Order Status: Accepted	46 Lower 36				Come From	
	45 Lower 35				梁永诊所	
	44 34					
	43 42 41 31 32 33					

Figure 22

In the sending interface, users can export case information to disk. Click the "Export" button, and the system will pop up a save dialog box, as shown in the figure. The user can set the export name, export path, and export format, and click "OK" to proceed.


Export		×
导出文件夹名称		
test		
Export Path		
D:\KP-Digital-Data\		
Export Data Format		
STL OBJ	PLY	
	Cancel	ОК

Figure 23

3.3.6 Open Case

In the case list, double-click the case to open the case details page.

3.3.7 Copy Case



In the case list, click the copy button to copy the specified case.

Figure 24

3.4 Patient Management Module

3.4.1 New Patient

Click the "New Patient" button and enter patient information in the pop-up dialog box. The name is a required field and the rest are optional fields. After entering, click the "OK" button to proceed.

			New Patient		×		50 🖗 🕹 🕸	ΘG
KP Digit	al	R)	Organization*	ID				
			梁永诊所	Please enter				
Case List Patient List	t Plugin List		Name*			G Patient Name	New Patie	nt
	-		Please enter			G		_
			Gender Man O Woman O Not	specified				
999	Male		Birthday	Age		2023-08-18 14:08:18		
999	Male		Please enter	Please enter		2023-08-18 14:08:05		
855	Male		Phone	Email		2023-08-18 12:02:57		
			Please enter	Please enter		2023-00-10 12:02:37		
888	Male		Address			2023-08-18 12:02:52		
777	Male		Please enter			2023-08-18 12:01:18		
π	Male		Medical history			2023-08-18 12:01:15		
666	Male		Please enter Remark			2023-08-18 11:51:45		
666	Male		Please enter			2023-08-18 11:51:45		
666	Male					2023-08-18 11:51:45		
666	Male			Cancel OK		2023-08-18 11:51:44		
							2 in total Show per page 20 = 🚺 1	/3

Figure 25

3.4.2 Search patients

In the patient list, users can enter the patient name for fuzzy search, and the search results are displayed in the list, as shown in the figure.

Case List Pa	tient List Plugin List			C Patient Nar	ne New Patient
Patient Name	Patient Gender	Patient Age	Contact Information	Creation time	Operation
999	Male			2023-08-18 14:08:18	Ø 📋
999	Male			2023-08-18 14:08:05	Ø 🗊
888	Male			2023-08-18 12:02:57	Ø 🗊
			Figure 26		



3.4.3 Modify Patient

In the patient list, click the edit button, enter the information that needs to be modified in the popup dialog box, and click "OK" to save the modified information, as shown in the figure:

		Modify Patient		×	
KP D	igital 🛛	Organization*	ID		
		梁永诊所	Please enter		
Case List Patie	ent List Plugin List	Name*		C Patient	New Patient
		999			
		Gender Man O Woman	O Not specified	Creation time	
999	Male	Birthday	Age	2023-08-18 14:08:18	
999	Male	Please enter	Please enter	2023-08-18 14:08:05	
855	Male	Phone Please enter	Email Please enter	2023-08-18 12:02:57	
888	Male	Address		2023-08-18 12:02:52	
777	Male	Please enter		2023-08-18 12:01:18	
m	Male	Medical history		2023-08-18 12:01:15	
666	Male	Remark		2023-08-18 11:51:45	
666	Male	Please enter		2023-08-18 11:51:45	
666	Male			2023-08-18 11:51:45	
666	Male		Cancel OK	2023-08-18 11:51:44	
					42 in total Show per page 20 = 🔇 1/3 🚺

Figure 27

3.4.4 Delete Patient

In the patient list, click the delete icon and a confirmation prompt will pop up. After the user confirms the deletion, the system will delete the selected patient information.




Figure 28

For patients who have already created a case, their information cannot be deleted.




Figure 29

3.5 Model Scan

In the case editing interface, click the "Descan" button at the top to enter the scanning interface, as shown in the figure:




Figure 30

3.5.1 Functional Area Description

1	Scanning steps: Display the current step during the scanning process
2	Upper and lower jaw/occlusal switch: switch to the current scanned jaw
3	3D model display area: displays the 3D model generated by scanning
4	Video preview area: Display real-time video of scanned teeth, and display scanner status when not scanning
5	Image toolbar: can select and manage 3D images
6	Soft tissue removal switch: When enabled, soft tissue can be automatically removed when scanning the mouth
7	High gloss surface switch: When scanning bright metal, this switch can be turned on for scanning
8	Mold scanning switch: When enabled, the mold area can be scanned. After exiting, the 3D model of the mold will be flipped and merged into the 3D model of the teeth
9	Option menu: allows for software configuration and viewing of software versions
10	Display toolbar: You can change the display of 3D models

3.5.2 Scan Icon Function Description

Scanning steps						
() Scan	Scan button for scanning the upper jaw, lower jaw, and occlusal relationship					
© Refine	Optimize button to refine the 3D model obtained from scanning and use various tools to check the optimization results					
Finish	Complete button, complete case information, and save scan results					
Scan section switching						
	Lower jaw button, switch to lower jaw for scanning or viewing					

	Upper jaw button, switch to upper jaw for scanning or viewing
	Bite button, switch to bite relationship for scanning or viewing
	Upper jaw scanning rod button, switch to upper jaw scanning rod scanning or viewing
	Lower jaw scanning rod button, switch to lower jaw scanning rod scanning or viewing
	Swap upper and lower jaw buttons to exchange scanned upper and lower jaw data
image tooll	par
»	Trim button, freely draw closed-loop curve to delete unnecessary data
3	Undoes the previous action
«	Return to the previous level of operation
6	Lock button to lock one or more areas on the model to prevent updates from subsequent scans
ଲ	Mark Teeth button to enter the Mark Teeth toolbar
8	Add tooth marker button
2	Delete Tooth Marker Button
S [℃]	Move tooth marker button
ದ್	Inverted check button to calculate the inverted teeth from the current perspective
€	Delete button to delete the generated model data
ă	Neck edge line positioning direction button, set the neck edge line positioning direction to the current viewing angle direction
\$	Neck edge line button, add neck edge line
\$	Neck edge line editing button
G	Neck edge line delete button
Show toolb	
	True color button, when selected, displays the 3D model in actual color, and when not selected, displays the 3D model in monochrome
0	Realistic display button, with a more realistic display effect when selected
5	Scan quality prompt button, displayed in green for areas with good scanning results
(Display the 3D center button. When selected, the rotation center of the current model will be displayed when rotating the 3D model
	Full screen display button to zoom the 3D model to the optimal size to fit the display range
6	Screen snapshot button to take a snapshot of the 3D model displayed on the screen
Ø	Oral endoscope button, select 2D intraoral image
Configurati	
00	Option menu button to open option settings and about dialog boxes

3.5.3 Option Configuration

						×
		\square	Ţ	\square	U	
-	General	Scan	Display	Save	Tools	
	Show instru	ictions				
			ок	Cancel		
Dicploy	nromnt. Whon	colocited the	Figur		l animations for	• some functions.
Display	prompt. when	selecteu, me	Soltware will	pop up tutoria		×
	(<u> </u>)			\sim	د ک	
	General	Scan	Display	Save	Tools	
	Hole highligh	t color		-		
		mote control		_		
	🥑 Enable sca	anning sound		Radar		•
	🥑 Enable wa	arning sound				

Scan position Occlusal pressure level OCCUSAL PRESSURE level OCCUSAL PRESSURE level OCCUSAL PRESSURE level OCCUSAL PRESSURE level Cancel

Figure 32

Hole marking color: The software fills the tooth holes in the model with the color specified during scanning and optimization


Enable remote sensing control: When selected, after scanning stops, the scanned 3D model can be rotated by rotating the scanner

Enable scanning sound: When selected, play a prompt sound when scanning images

Enable warning sound: When selected, if the scanning time of the image exceeds the recommended threshold, or if there is strong light during scanning, or if the scanning performance decreases, a warning sound will be played

Scanning pose: switching the direction displayed during model scanning

Degree of occlusal pressure: Adjust the occlusal pressure





Interface element size: Select the size of the interface icon to adapt to different screen resolutions Monochrome mode color: Select the color when displaying a 3D model in monochrome

General	Scan	D isplay	C Save	O Tools	
< Save bite c	lata				
ID orthod	ontic data				
O HD color					
		ОК	Cancel		

 \times

 \times



Save bite data: When selected, save the bite file

High definition orthodontic data: When selected, the orthodontic scan file will be saved in high resolution

High definition color: When selected, the map is saved in high-definition format

	\square	Ţ		ନ	
General	Scan	Display	Save	Tools	
Shut down the	e wired scanner a	fter	10 minutes		•
Shut down the	e wireless scanner	rafter	3 minutes		•
		Create sup	oport file		
		ок	Cancel		

Figure 35

Turn off wired scanner after this time: Users can customize the automatic shutdown time of the scanner Creating supporting files: exporting software and firmware logs

3.6 Model Preview



Figure 36

3.6.1 Functional Area Description

1	Model file list: All scan generated model files are listed here
2	Angle adjustment bar: adjust the display direction of the model
3	Model display area
4	Show toolbar

3.6.2 Model Preview Icon Function Description

0	Show/Hide Model Button
Û	Delete button to delete the current model data
	Front view button, adjust the current perspective of the model to the front view direction
	Back view button, adjust the current perspective of the model to the back view direction
And the second second	Top view button, adjust the current perspective of the model to the top view direction
5999 B	Look up button, adjust the current perspective of the model to look up direction
B	Right view button, adjust the current perspective of the model to the right view direction
BB	Left view button, adjust the current perspective of the model to the left view direction
۲	Gypsum color button, switch model color to gypsum color
	Angle adjustment button to adjust the current model size to the optimal display size

3.7 Plug-in Manager

3.7.1 Search plugin

Enter the plugin name, click the search button, and query for the specified plugin

Case List Patient List	Plugin List		C Smile
SmileScan V1.0.0 2023-07-26			
	Pre-Installed		

Figure 37

3.7.2 View plugin information

Click on the plugin to view the plugin information

KP Digital	SmileScan	KP Digital		🛜 ⋤ 💮 靠 liangyong 🕞 🕞
Case List Patient List Plugin List	V1.0.0 2023-07-26	P	Pre-Installed	C Plugin Name
SmileScan V100 2023-07-76	Description SmileScan is a scanning plugin			
Pre-in				
	Version	Required KP Digital Version		
	V1.0.0	V1.0.0	_	
	Initial Release Date	Latest Release Date		
	2023-07-26	2023-07-26		
			Close	

Figure 38

3.7.3 Install plug-in

Click on the plugin and click the install button on the plugin information page to install it

3.7.4 Uninstalling plugins

Click on the installed plugin and click the uninstall button on the plugin information page to uninstall it

3.8 System settings

Click on the settings icon in the upper right corner of the software to enter the settings interface, which includes common settings, help, and software update functions.

3.8.1 Common

Common settings include functions such as multilingual switching, exporting folder name format settings, storage path settings, backup path settings, and export path settings.

Setting		×
Common	Language	
Help	English Export folder name format	
Update Software	Patient Name	
	Storage Path D:\KP-Digital-Data\	
	Backup Path D:\KP-Digital-Data\backup\	
	Export Path	
	D:\KP-Digital-Data\	ОК

Figure 39

3.8.2 Help

The help interface displays the name of the software development unit and company website. Click the "Open" button to display the software user manual.

Setting		×
Common	Manufacturer: Guilin Kevin Peter Technology Co., Ltd.	
Help	Company website: http://www.kevinpetertech.com	
Update Software	Software manual:	

Figure 40


3.8.3 Update software

The software update interface displays the version information of the current software. Clicking on "Check for updates" will connect to the cloud to check for new versions.

×

Figure 41

When the current version is the latest:





When there is a new version, you will be prompted whether to download the upgrade package. Click "Yes" to start downloading the upgrade package. After downloading, you will be prompted whether to install it. Click "Yes" to immediately install it.




Figure 43




Figure 44

3.9 Connection Assistant

The mobile app needs to connect to the server to view case data. The mobile network and computer network are in the same local area network. Click on the connection assistant, and the mobile app can scan the QR code on the interface to log in. You can also manually enter the IP and Port on the connection assistant to log in.

Connection Assistant					×	
Scanning connnection on mobile app	М	anual co	nnection	on mobi	le app	
MIKE	IP	10	66	66	226	
No. 10		172	20	20	226	
23,215,2	Port		80	30		
Scan QR code for connection on the KP Digital mobile app			s and port to connect	number in	the KP	
				R	efresh	

Figure 45

3.10 User Management

After logging into the software, you can click the My Information button to view the basic information

of the logged in user, edit user information, or click the Exit Login button to return to the login page.

My Information	×
Organization Name	
梁永诊所	
Email Address	
657587346@qq.com	
Name 梁永	
Edit Log Out Close	
Figure 46	

3.11 Cloud management

Click on the "KP Digital Cloud Platform" button to go to the cloud platform webpage for cloud management

KF V1.0.0	P Digital	New Case	KP Digital Cloud Platform	~	Fo 😳 😩 🗫	O C
Case List	Patient List Plugin List				Case No. , Patient Name	8
Case No.	Patient Name	Case Type	Case Status	Last Modified Time	Operation	
006	999	Repair	Waiting to scan	2023-08-19 16:08:45	🖻 📋	
0041	hhhhhh慢嘎嘎嘎嘎嘎!11	Repair	Waiting to deliver	2023-08-19 07:16:21		
006	999	Repair	Waiting to scan	2023-08-18 14:17:58	🕒 🤠	
005	55	Repair	Waiting to scan	2023-08-18 11:41:10	🗊 🗊	
004	34	Repair	Waiting to scan	2023-08-18 11:16:48	()	
002	test	Repair	Accepted	2023-08-17 14:52:49	ڭ 🖸 🖆	
001	liangy	Repair	Accepted	2023-08-17 14:23:04		

7 in total Show per page 20 = (1/1)

Figure 47

4 Scope of application

As a software component of the oral digital impression instrument, its main function is to receive image information from the host for generating 3D models and managing patient orders.

This device is only allowed for professional and trained personnel to use. Proper use of this device will not cause side effects; If this device is not used correctly, it will result in incorrect image acquisition.