

**Conversion Chart ParaPost & FiberKleer > Mecodont Post**

ParaPost	ParaPost Plus	ParaPost XP	Color	Ø inches	Ø mm	Post #	FiberKleer	Color	Ø mm	Mecodont	Color	Ø inches	Ø mm	Post #	
<b>FIBER POST INTRO.KITS</b>															
PF160	PF170	Glass Fiber parallel sided & serrated, including: 10 posts ( 3 x 1.14, 3 x 1.25, 2 x 1.40, 2 x 1.5mm), 4 drills (1 of each size), 1 universal driver					N83A	FibreKleer 4X Parallel Posts5 ea. of sizes 1.00mm red, 1.25mm yellow, 1.50mm blue, 3 Drills- 1 ea. of sizes			09-100	Glass Fiber parallel sided & serrated including: 21 posts (3 of each size) 7 drills (1 of each size), 1 vial of endo stops			
<b>FIBER POST REFILLS</b>															
PF1613	PF1713		Brown	0.36	0.90	3	Size Not Available			09-101	Orange	0.36	0.90	3	
PF1614	PF1714		Yellow	0.40	1.00	4	N83AA	Red	1.00	09-102	Pink	0.40	1.00	4	
PF16145	PF17145		Blue	0.45	1.14	4.5	N83CB	Blue	1.14	09-103	Grey	0.45	1.14	4.5	
PF1615	PF1715		Red	0.50	1.25	5	N83AB	Yellow	1.25	09-104	Beige	0.50	1.25	5	
PF16 55	PF1755		Purple	0.55	1.40	5.5	N83BB	Purple	1.38	09-105	Purple	0.55	1.40	5.5	
PF1616	PF1716		Black	0.60	1.50	6	N83AC	Blue	1.50	09-106	White	0.60	1.50	6	
ParaPost does not have equivalent size to Mecodont							Size Not Available			09-107	Black	0.70	1.75	7	
<b>TITANIUM POST INTRO.KITS</b>															
		P780T	Titanium Parallel sided & vented, including: 25 posts, 1 of each size, 1 vial of endo stops				FiberKleer not available in Titanium			09-100	Titanium Parallel sided & vented, including: 21 posts (3 of each size) 7 drills (1 of each size), 1 vial of endo stops				
<b>TITANIUM POSTS REFILLS</b>															
P8843	P2843	P7843	Brown	0.36	0.90	3				09-101	Orange	0.36	0.90	3	
P844	P2844	P7844	Yellow	0.40	1.00	4				09-102	Pink	0.40	1.00	4	
P8445	P2845	P78445	Blue	0.45	1.14	4.5				09-103	Grey	0.45	1.14	4.5	
P845	P2845	P7845	Red	0.50	1.25	5				09-104	Beige	0.50	1.25	5	
P8455	P28455	P78455	Purple	0.55	1.40	5.5				09-105	Purple	0.55	1.40	5.5	
P846	P2846	P7846	Black	0.60	1.50	6				09-106	White	0.60	1.50	6	
P847	P2847	P7847	Green	0.70	1.75	7				09-107	Black	0.70	1.75	7	
<b>STAINLESS STEEL POST INTRO.KITS</b>															
		P780	Stainless Steel Parallel sided & vented, including: 25 posts, 1 of each size, 1 vial of endo stops				FiberKleer not available in Stainless Steel			09-001	Stainless Steel Parallel sided & vented, including: 21 posts (3 of each size) 7 drills (1 of each size), 1 vial of endo stops				
<b>STAINLESS STEEL POST REFILLS</b>															
P443	P2443	P7443	Brown	0.36	0.90	3				09-002	Orange	0.36	0.90	3	
P444	P2444	P7444	Yellow	0.40	1.00	4				09-003	Pink	0.40	1.00	4	
P4445	P24445	P74445	Blue	0.45	1.14	4.5				09-004	Grey	0.45	1.14	4.5	
P445	P2445	P7445	Red	0.50	1.25	5				09-005	Beige	0.50	1.25	5	
P4455	P24455	P74455	Purple	0.55	1.40	5.5				09-006	Purple	0.55	1.40	5.5	
P446	P2446	P7446	Black	0.60	1.50	6				09-007	White	0.60	1.50	6	
P447	P2447	P7447	Green	0.70	1.75	7				09-008	Black	0.70	1.75	7	
<b>DRILLS REFILLS</b>															
P423	P6230		Brown	0.36	0.90	3	Size Not Available			09-200	Orange	0.36	0.90	3	
P424	P6240		Yellow	0.40	1.00	4	N18DA	Red	1.00	09-201	Pink	0.40	1.00	4	
P4245	P6245		Blue	0.45	1.14	4.5	N18DD	Purple	1.125	09-202	Grey	0.45	1.14	4.5	
P425	P6250		Red	0.50	1.25	5	N18DB	Yellow	1.25	09-203	Beige	0.50	1.25	5	
P42-55	P6255		Purple	0.55	1.40	5.5	N18DE	Green	1.38	09-204	Purple	0.55	1.40	5.5	
P426	P6260		Black	0.60	1.50	6	N18DC	Blue	1.50	09-205	White	0.60	1.50	6	
P427	P6270		Green	0.70	1.75	7	Size Not Available			09-206	Black	0.70	1.75	7	