



## Main Differences Between Woodpecker PT-B and PT-A Scaler & Air Polisher Combo Models

| No. | Feature                   | PT-B  | PT-A  |
|-----|---------------------------|---|---|
| 1   | Appearance                |  |  |
| 2   | Screen Parameters         | Higher Resolution Screen<br>1280*720"   | Lower Resolution Screen<br>1024*600   |
| 3   | Heating Power & Functions | 300W Heating Power<br>Allows for Auto-Preheating<br>4 Temperature Settings        | 50W Heating Power<br>No Preheating Function<br>No Heating Level Options             |
| 4   | Handpiece                 | Magnetic induction type   | Mechanical button type  |
| 5   | Foot Pedal                | Wireless multi-function Foot Pedal  | Wired Multi-Function Foot Pedal   |
| 6   | Power Supply              | Built-In Switch-Mode Power Supply   | External Power Transformer  |
| 7   | Pricing                   | \$6,497.00  | \$3,997.00  |